



800-872-1001

Accession No. [REDACTED]  
 Received 11/18/2015  
 Reported 11/18/2015 03:40 AM

Doctor [REDACTED]

Dz: Liver  
 Wk: Heart  
 ^CA

Owner **ESNES** Pet Name **SHELBY** Species **Canine** Breed **Mixed** Sex **SF** Pet Age **12Y** Chart# **N**

Superchem				Complete Blood Count			
Tests	Results	Ref. Range	Units	Tests	Results	Ref. Range	Units
Total Protein	6.4	5.0-7.4	g/dL	WBC	7.6	4.0-15.5	10 <sup>3</sup> /μL
Albumin	4.0	2.7-4.4	g/dL	RBC	6.5	4.8-9.3	10 <sup>6</sup> /μL
Globulin	2.4	1.6-3.6	g/dL	Hemoglobin	15.3	12.1-20.3	g/dL
A/G Ratio	1.7	0.8-2.0	Ratio	Hematocrit	47	36-60	%
AST (SGOT)	29	15-66	U/L	MCV	72	58-79	fL
<b>M</b> ALT (SGPT)	<b>153 (HIGH)</b>	<b>12-118</b>	<b>U/L</b>	MCH	23.5	19-28	pg
<b>A</b> Alk Phosphatase	124	5-131	U/L	MCHC	32.6	30-38	g/dL
<b>A</b> GGTP	5	1-12	U/L	Platelet Count	<b>A</b> 371	170-400	10 <sup>3</sup> /μL
Total Bilirubin	0.1	0.1-0.3	mg/dL	Platelet EST	Adequate	Adequate	
<b>A</b> Urea Nitrogen	20	6-31	mg/dL	<b>Differential</b>	<b>Absolute</b>	<b>%</b>	
Creatinine	1.3	0.5-1.6	mg/dL	Neutrophils	5092	67	2060-10600 /uL
BUN/Creatinine Ratio	15	4-27	Ratio	Bands	0	0	0-300 /uL
Phosphorus	4.1	2.5-6.0	mg/dL	Lymphocytes	1900	25	690-4500 /uL
Glucose	93	70-138	mg/dL	Monocytes	304	4	0-840 /uL
<b>A</b> Calcium	11.4	8.9-11.4	mg/dL	Eosinophils	304	4	0-1200 /uL
Magnesium	1.6	1.5-2.5	mEq/L	Basophils	0	0	0-150 /uL
Sodium	146	139-154	mEq/L				
Potassium	4.8	3.6-5.5	mEq/L				
Na/K Ratio	30						
Chloride	111	102-120	mEq/L				
Cholesterol	261	92-324	mg/dL				
Triglycerides	107	29-291	mg/dL				
Amylase	523	290-1125	U/L				
PrecisionPSL™	80	24-140	U/L				

Pancreatitis is unlikely, but a normal PrecisionPSL™ result does not completely exclude pancreatitis as a cause for gastrointestinal signs.

**New from Antech Diagnostics!!!!**

**Now you never have to order a PLI again! Effective immediately Antech Diagnostics introduces the Precision Pancreatic Sensitive Lipase (PrecisionPSL™) assay replacing the traditional lipase assay at NO additional cost!**

Antech's new PrecisionPSL™ test is replacing the traditional lipase assay on all Superchem profiles. PrecisionPSL™ is analogous in performance to pancreatic lipase immunoreactivity concentrations, precluding the need to perform PLI-type testing. Interpretive comments will accompany each report. There is no additional charge for this new lipase test. The previous lipase assay is no longer available.

<b>V</b> CPK	63	59-895	U/L
Comment	Hemolysis 1+ No significant interference.		

Test Requested	Results	Reference Range	Units
<b>V</b> LDH			
LDH	72	20-500	U/L
Comment	Hemolysis 1+ No significant interference.		